


# Real Time PCR

TOP QUALITY MOLECULAR DIAGNOSTICS  
FOR INFECTIOUS DISEASES



All Products Manufactured & Developed by Dia.Pro Italy

*Your IVD  
specialists!*



## HEPATITIS

KIT NAME	DESCRIPTION
<b>Hepatitis B</b>	
HBV DNA (QT)	Quantitative detection of Hepatitis B Virus DNA
<b>Hepatitis Delta</b>	
HDV RNA	Quantitative detection of Hepatitis D Virus RNA

**KIT FORMAT: 25, 50, 100, 150 reactions**

## TRANSPLANTATION

KIT NAME	DESCRIPTION
<b>Cytomegalovirus</b>	
CMV DNA (QT) 2 <sup>nd</sup> Gen.	Quantitative detection of Cytomegalovirus DNA
<b>Herpes Simplex Virus 1</b>	
HSV1 DNA (QT)	Quantitative of Human herpes simplex virus 1 DNA
<b>Herpes Simplex Virus 2</b>	
HSV2 DNA (QT)	Quantitative of Human herpes simplex virus 2 DNA
<b>Herpes Virus 6</b>	
HHV6 DNA (QT)	Quantitative detection of Herpesvirus 6 DNA
<b>Herpes Virus 8</b>	
HHV8 DNA (QT)	Quantitative detection of Herpesvirus 8 DNA
<b>Epstein Barr Virus</b>	
EBV DNA (QT)	Quantitative detection of Epstein-Barr Virus DNA
<b>Varicella Zoster Virus</b>	
VZV DNA (QT)	Quantitative detection of Varicella-Zoster Virus DNA
<b>BK Virus</b>	
BKV DNA (QT)	Quantitative detection of BK Virus DNA
<b>JC Virus</b>	
JCV DNA (QT)	Quantitative detection of JC Virus DNA
<b>Enterovirus spp</b>	
ENTEROVIRUS RNA	Qualitative detection of Enterovirus RNA

**KIT FORMAT: 25, 50, 100, 150 reactions**

## MISCELLANY

KIT NAME	DESCRIPTION
<b>Toxoplasma Gondii</b>	
TOXOPLASMA GONDII DNA	Qualitative detection of Toxoplasma Gondii DNA
<b>Legionella Pneumophila</b>	
LEGIONELLA PENUMOPHILA DNA	Qualitative detection of Legionella Pneumophila DNA
<b>Trypanosoma Cruzi (Chagas Disease)</b>	
TRYPANOSOMA CRUZI DNA	Qualitative detection of Trypanosoma cruzi DNA
TRYPANOSOMA CRUZI DNA (QT)	Quantitative detection of Trypanosoma cruzi DNA
<b>Chlamydia trachomatis</b>	
CHLAMYDIA TRACHOMATIS DNA	Qualitative detection of Chlamydia trachomatis DNA

**KIT FORMAT: 25, 50, 100, 150 reactions**

## ACCESSORY KIT

KIT NAME	DESCRIPTION
<b>RNA Retrotranscription Reagents</b>	
RNA Retrotranscription Kit	Reverse transcription to cDNA of highly purified RNA coming from virus (vRNA)
HDV RNA Retrotranscription Kit	Reverse transcription to cDNA of Hepatitis delta virus RNA

**KIT FORMAT: 25, 50, 100, 150 reactions**

## ASSAY PROTOCOLS VALIDATED ON

- **ABI 7500 Sequence Detection System** (Software SDS version 1.3.1, Applied Biosystems™)
- **CFX 96** (Software CFX Manager version 1.7 Biorad™)
- **MX3000P** (Software MxPro version 4.01, Stratagene/Agilent Technologies)
- **ROTO-GENE Q Series** (Qiagen™)

**Assay protocols may be adapted on other RT PCR platforms upon request**

# TOP SENSITIVITY AND SPECIFICITY

Your IVD  
specialists!

## MAIN FEATURES

- ✓ CE certified
- ✓ Kits contains all reagents necessary to carry out test
- ✓ Internal control to verify DNA extraction and amplification
- ✓ Quantitative tests: standard curve on 4 points
- ✓ High sensitivity and specificity
- ✓ Small Reaction volume of 25  $\mu$ l: 20  $\mu$ l Reaction Mix + 5  $\mu$ l samples
- ✓ Sample: Blood, Serum and Plasma, Liquor, Biopsies
- ✓ Storage temperature: +2/8°C (no shipment in dry ice needed)
- ✓ Shelf life: 12 months at +2/8°C
- ✓ Kit format: 25, 50, 100, 150 reactions: non waste

**MADE IN ITALY**

For information visit: [www.diapro.it](http://www.diapro.it)



DIA.PRO Diagnostic Bioprobes S.r.l.  
Via Carducci, 27 - 20099 Sesto San Giovanni (MI) - Italy  
Tel. +39 02 27.00.71.61/64.50 Fax +39 02 44.38.67.71  
<http://www.diapro.it> - email: [info@diapro.it](mailto:info@diapro.it)

**EN ISO 13485:2013 Certified**