



DIA. PRO

PRODUCT
CATALOGUE
2022



Your IVD
Specialists!



CLIA

PCR

ELISA

www.diapro.it





Made in Italy, distributed worldwide.

Dia.Pro Diagnostic Bioprobes srl is a Global key player in *in vitro* diagnostics for 25 years.

Dia.Pro manufactures and markets a wide range of high-quality infectious disease assays and automated solutions worldwide.

Our commitment is to enable our customers around the world provide efficient and reliable diagnosis for patients based on clinical value.

All our in vitro products (Reagents and Automated solutions) are CE marked and strictly comply with requirements of the European in vitro Diagnostic Directive (IVDD) 98/79/EC.

The complete manufacturing process of Dia.Pro's products, from the Production of Raw Material to the Quality Control Release, is carried out in our 13485 certified facility (3,000 square meters) located in Milan, Italy.

Our deliveries are managed efficiently, as to answer our customers' needs, by supplying complete orders. We are recognised for our flexible and competitive prices, as well as our commitment to excellence in Customer Service.



SUMMARY



CE

CLIA

CLIA SYSTEM	7
TECHNICAL SPECIFICATIONS	8
TEST MENU	9
COVID-19	9
BLOOD SCREENING	10
TORCH	10
AUTOIMMUNITY	11
SEROLOGY	11
EMERGING TROPICAL DISEASES	11

SARA



Bench-top
Fully Automated
RANDOM ACCESS SYSTEM
Based on Magnetic Particles

TECHNICAL SPECIFICATIONS

Description	A fully automated bench top system with continuous random access that offers an expanding menu of over 50 immunoassays that can be run side by side
Throughput	Max. 120 test per hour in random access
Time to first Result	15 minutes (assay dependent)
Walkaway Capacity	400 test; auto dilution; auto repeat-testing; auto reflex testing
Dispensing System	Zero carryover with disposable Tip System; sample integrity checking, including low volume sample and clot detection
Methodologies Used	Chemiluminescence Immunoassays (CLIA)
Startup Time	5 minutes startup: 15 minutes daily maintenance with no hands on time
Sample Handling	On board 90 sample racks with continuous feed capability; Barcode identification; primary, secondary and pediatric tubes: sample volume 10-200µl
Reagents Onboard	16 Refrigerated reagents with barcode Identification RFID; 50-200 tests/kit (assay dependent); extensive 4 weeks on board stability; automatic inventory, residency tracking and flagging; calibration frequency of 4 weeks validity tracked and flagged

SOFTWARE SPECIFICATION

Software Platform	Minimized operator intervention using Window-based fully multitasking software intuitive multilanguage interface for easy navigation; touch screen monitor
Quality Control	QC timeline using colors and symbols to indicate changes in lot numbers; user-definable Westgard rules and Z-scores for comprehensive QC management; integrated data management; full patients search by ID, name, data and test; session and patient report printing
Host Connection	Bidirectional to any LIS in conformance with ASTM E 1394 and ASTM EI 381; networks via Ethernet (10/100 Base T) or WiFi.

INSTRUMENT SPECIFICATION

Dimension and Weight	Height 90 cm - Width 190 cm - Depth 78 cm - Weight +/- 200 kg
Operating Conditions	Temperature 15-32°C; power supply 100-120 (200-240) VAC; 50-60 Hz; consumption 600 VA
Certification	ISO, EEC, IVD, CE; detailed specifications available on request

COVID-19



CODE	KIT NAME	DESCRIPTION	TEST
RACOV19M.CE	DIA.CHEMILUX SARS-CoV-2NCP IgM	Determination of IgM antibodies	100
RACOV19A.CE	DIA.CHEMILUX SARS-CoV-2NCP IgA	Determination of IgA antibodies	100
RACOV19G.CE	DIA.CHEMILUX SARS-CoV-2NCP IgG	Determination of IgG antibodies	100
RACOV19SPIKEAB.CE	DIA.CHEMILUX SARS-CoV-Spike/RBD Ab	Determination of IgG antibodies to COVID-19 Spike/RBD antigens	100
RACOV19NEUTRAB.CE	SARS-CoV-2 ACE2-RBD NEUTRALIZATION	Determination of the neutralization activity of antibodies to SARS-CoV-2, induced by natural infection and/or vaccination, able to block the interaction of Receptor Binding Domain (RBD) and ACE2	100



BLOOD SCREENING

CODE	KIT NAME	DESCRIPTION	TEST	NOTES
RASAG.CE	DIA.CHEMILUX HBs Ag		100	
RAIVCOMB.CE	DIA.CHEMILUX HIV Ab&Ag	Qualitative	100	
RACVAB.CE	DIA.CHEMILUX HCV Ab	Qualitative	100	
RASIAB.CE	DIA.CHEMILUX Syphilis Ab	Qualitative	100	
RATCAB.CE	DIA.CHEMILUX T cruzi Ab	Qualitative	100	
RAHTLVAB.CE	DIA.CHEMILUX HTLV I&II Ab	Qualitative	100	
RABCAB.CE	DIA.CHEMILUX HBC Total Ab	Quantitative	100	



TORCH

CODE	KIT NAME	DESCRIPTION	TEST	NOTES
RATOXOG.CE	DIA.CHEMILUX Toxoplasma IgG	Quantitative	100	
RATOXOM.CE	DIA.CHEMILUX Toxoplasma IgM	Quantitative	100	IgM Capture System
RARUBG.CE	DIA.CHEMILUX Rubella IgG	Quantitative	100	
RARUBM.CE	DIA.CHEMILUX Rubella IgM	Quantitative	100	IgM Capture System
RACMVG.CE	DIA.CHEMILUX Cytomegalovirus IgG	Quantitative	100	
RACMVM.CE	DIA.CHEMILUX Cytomegalovirus IgM	Quantitative	100	IgM Capture System
RAHSV1G.CE	DIA.CHEMILUX HSV1 IgG	Quantitative	100	
RAHSV1M.CE	DIA.CHEMILUX HSV1 IgM	Qualitative	100	IgM Capture System
RAHSV2G.CE	DIA.CHEMILUX HSV2 IgG	Quantitative	100	
RAHSV2M.CE	DIA.CHEMILUX HSV2 IgM	Qualitative	100	IgM Capture System
RAHSVG.CE	DIA.CHEMILUX HSV1&2 IgG	Quantitative	100	
RAHSVM.CE	DIA.CHEMILUX HSV1&2 IgM	Qualitative	100	IgM Capture System

AUTOIMMUNITY



CODE	KIT NAME	DESCRIPTION	TEST	NOTES
RASM.CE	DIA.CHEMILUX IgG anti-sm		100	
RASSA52.CE	DIA.CHEMILUX IgG anti-SSA52	Quantitative	100	
RASSA60.CE	DIA.CHEMILUX IgG anti-SSA60	Quantitative	100	
RASSB.CE	DIA.CHEMILUX IgG anti-SSB	Quantitative	100	
RARNP.CE	DIA.CHEMILUX IgG anti-RNP	Quantitative	100	
RAJO1.CE	DIA.CHEMILUX IgG anti:JO1	Quantitative	100	
RASCL70.CE	DIA.CHEMILUX IgG anti-SCL70	Quantitative	100	

SEROLOGY



CODE	KIT NAME	DESCRIPTION	TEST	NOTES
RATETG.CE	DIA.CHEMILUX TOXOID IgG	Quantitative	100	
RAEVAB.CE	DIA.CHEMILUX HEV Ab	Qualitative	100	
RADAB.CE	HEPATITIS D VIRUS	Qualitative	100	
RAMALAB.CE	MALARIA	Quantitative	100	

EMERGING TROPICAL DISEASES



CODE	KIT NAME	DESCRIPTION	TEST	NOTES
RAZIKVAB.CE	DIA.CHEMILUX ZIKV Ab	Qualitative	100	
RACHIKVAB.CE	DIA.CHEMILUX CHIKV Ab	Qualitative	100	
RADENGVAB.CE	DIA.CHEMILUX Dengue Virus Ab	Qualitative	100 Test/kit	

COMING SOON

HAV Ab • HBs Ab • HDV IgM • dsDNA Ab • Transglutaminase IgA&IgG • Gliadin deamidated IgA&IgG

REAL TIME PCR

TEST MENU

COVID-19	13
HEPATITIS	13
TRANSPLANTATION	14
MISCELLANY	17
ACCESSORY KIT	18
PLATFORM VALIDATION	18

COVID-19

COVID-19

SARS-CoV-2

CODE	KIT NAME	DESCRIPTION	TEST
COV19RNALYO.CE	COVID-19 RNA Vs 2	Lyophilized One Step Multiplex Real Time RT PCR mix for the detection of novel Coronavirus SARS-CoV-2 RNA in human respiratory tract specimens	100
COV19RNALYO.CE.200	COVID-19 RNA Vs 2	Lyophilized One Step Multiplex Real Time RT PCR mix for the detection of novel Coronavirus SARS-CoV-2 RNA in human respiratory tract specimens	200

HEPATITIS



HEPATITIS B

CODE	KIT NAME	DESCRIPTION	TEST
HBVDNAQT.CE.25	HBV DNA	Quantitation (QT)	25
HBVDNAQT.CE	HBV DNA	Quantitation (QT)	50
HBVDNAQT.CE.100	HBV DNA	Quantitation (QT)	100
HBVDNAQT.CE.150	HBV DNA	Quantitation (QT)	150

HEPATITIS DELTA

CODE	KIT NAME	DESCRIPTION	TEST
DRNA.CE.25	HDV RNA	Quantitation (QT)	25
DRNA.CE	HDV RNA	Quantitation (QT)	50
DRNA.CE.100	HDV RNA	Quantitation (QT)	100
DRNA.CE.150	HDV RNA	Quantitation (QT)	150



TRANSPLANTATION

CYTOMEGALOVIRUS

CODE	KIT NAME	DESCRIPTION	TEST
CMVDNAQT.2G.CE.25	CMV DNA	Quantitotition (QT) 2nd Generation	25
CMVDNAQT.2G.CE	CMV DNA	Quantitotition (QT) 2nd Generation	50
CMVDNAQT.2G.CE.100	CMV DNA	Quantitotition (QT) 2nd Generation	100
CMVDNAQT.2G.CE.150	CMV DNA	Quantitotition (QT) 2nd Generation	150

HERPES SIMPLEX VIRUS 1

CODE	KIT NAME	DESCRIPTION	TEST
HSV1DNAQT.CE.25	HSV1 DNA	Quantitation (QT)	25
HSV1DNAQT.CE	HSV1 DNA	Quantitation (QT)	50
HSV1DNAQT.CE.100	HSV1 DNA	Quantitation (QT)	100
HSV1DNAQT.CE.150	HSV1 DNA	Quantitation (QT)	150

HERPES SIMPLEX VIRUS 2

CODE	KIT NAME	DESCRIPTION	TEST
HSV2DNAQT.CE.25	HSV2 DNA	Quantitation (QT)	25
HSV2DNAQT.CE	HSV2 DNA	Quantitation (QT)	50
HSV2DNAQT.CE.100	HSV2 DNA	Quantitation (QT)	100
HSV2DNAQT.CE.150	HSV2 DNA	Quantitation (QT)	150

HERPES VIRUS 6

CODE	KIT NAME	DESCRIPTION	TEST
HHV6DNAQT.CE.25	HHV6 DNA	Quantitation (QT)	25
HHV6DNAQT.CE	HHV6 DNA	Quantitation (QT)	50
HHV6DNAQT.CE.100	HHV6 DNA	Quantitation (QT)	100
HHV6DNAQT.CE.150	HHV6 DNA	Quantitation (QT)	150

TRANSPLANTATION



HERPES VIRUS 8

CODE	KIT NAME	DESCRIPTION	TEST
HHV8DNAQT.CE.25	HHV8 DNA	Quantitation (QT)	25
HHV8DNAQT.CE	HHV8 DNA	Quantitation (QT)	50
HHV8DNAQT.CE.100	HHV8 DNA	Quantitation (QT)	100
HHV8DNAQT.CE.150	HHV8 DNA	Quantitation (QT)	150

EPSTEIN BARR VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
EBVDNAQT.CE.25	EBV DNA	Quantitation (QT)	25
EBVDNAQT.CE	EBV DNA	Quantitation (QT)	50
EBVDNAQT.CE.100	EBV DNA	Quantitation (QT)	100
EBVDNAQT.CE.150	EBV DNA	Quantitation (QT)	150

VARICELLA ZOSTER VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
VZVDNAQT.CE.25	VZV DNA	Quantitation (QT)	25
VZVDNAQT.CE	VZV DNA	Quantitation (QT)	50
VZVDNAQT.CE.100	VZV DNA	Quantitation (QT)	100
VZVDNAQT.CE.150	VZV DNA	Quantitation (QT)	150

BK VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
BKVDNAQT.CE.25	BKV DNA	Quantitation (QT)	25
BKVDNAQT.CE	BKV DNA	Quantitation (QT)	50
BKVDNAQT.CE.100	BKV DNA	Quantitation (QT)	100
BKVDNAQT.CE.150	BKV DNA	Quantitation (QT)	150



TRANSPLANTATION

JC VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
JCVDNAQT.CE.25	JCV DNA	Quantitotition (QT) 2nd Generation	25
JCVDNAQT.CE	JCV DNA	Quantitotition (QT) 2nd Generation	50
JCVDNAQT.CE.100	JCV DNA	Quantitotition (QT) 2nd Generation	100
JCVDNAQT.CE.150	JCV DNA	Quantitotition (QT) 2nd Generation	150

ENTEROVIRUS SPP

CODE	KIT NAME	DESCRIPTION	TEST
ENTERORNA.CE.25	ENTEROVIRUS RNA		25
ENTERORNA.CE	ENTEROVIRUS RNA		50
ENTERORNA.CE.100	ENTEROVIRUS RNA		100
ENTERORNA.CE.150	ENTEROVIRUS RNA		150

MISCELLANY

TOXOPLASMA GONDII

CODE	KIT NAME	DESCRIPTION	TEST
TOXDNA.CE.25	Toxoplasma Gondii DNA		25
TOXDNA.CE	Toxoplasma Gondii DNA		50
TOXDNA.CE.100	Toxoplasma Gondii DNA		100
TOXDNA.CE.150	Toxoplasma Gondii DNA		150

LEGIONELLA PNEUMOPHILA

CODE	KIT NAME	DESCRIPTION	TEST
LEPDNA.CE.25	LEGIONELIA PNEUMOPHILA DNA		25
LEPDNA.CE	LEGIONELIA PNEUMOPHILA DNA		50
LEPDNA.CE.100	LEGIONELIA PNEUMOPHILA DNA		100
LEPDNA.CE.150	LEGIONELIA PNEUMOPHILA DNA		150

TRYPANOSOMA CRUZI (CHAGAS DISEASE)

CODE	KIT NAME	DESCRIPTION	TEST
TCRUZIDNA.CE.25	TRYPANOSOMA CRUZI DNA		25
TCRUZIDNA.CE	TRYPANOSOMA CRUZI DNA		50
TCRUZIDNA.CE.100	TRYPANOSOMA CRUZI DNA		100
TCRUZIDNA.CE.150	TRYPANOSOMA CRUZI DNA		150
TCRUZIDNAQT. CE.25	TRYPANOSOMA CRUZI DNA	Quantitation (QT)	25
TCRUZIDNAQT. CE	TRYPANOSOMA CRUZI DNA	Quantitation (QT)	50
TCRUZIDNAQT. CE.100	TRYPANOSOMA CRUZI DNA	Quantitation (QT)	100
TCRUZIDNAQT. CE.150	TRYPANOSOMA CRUZI DNA	Quantitation (QT)	150

CHLAMYDIA TRACHOMATIS

CODE	KIT NAME	DESCRIPTION	TEST
CTDNA.CE.25	Chlamydia Trachomatis DNA		25
CTDNA.CE	Chlamydia Trachomatis DNA		50
CTDNA.CE.100	Chlamydia Trachomatis DNA		100
CTDNA.CE.150	Chlamydia Trachomatis DNA		150

ACCESSORY KIT

RNA RETROTRANSCRIPTION REAGENTS

CODE	KIT NAME	DESCRIPTION	TEST
RNART.CE.25	RNA	Retrotranscription KIT	25
RNART.CE	RNA	Retrotranscription KIT	50
RNART.HDV.CE.25	HDV RNA	Retrotranscription KIT	25
RNART.HDV.CE	HDV RNA	Retrotranscription KIT	50
RNART.HDV.CE.100	HDV RNA	Retrotranscription KIT	100

Assay protocols validated on



ABI 7500 Sequence Detection System

(Software SDS version 1.3. I, Applied Biosystems™)

CFX 96

(Software CFX Manager version I .7 Biorad™)

MX3000P

(Software MxPro version 4.0 I, Stratagene/ Agilent Technologies)

ROTOR-GENE Q Series

(Qiagen™)

Assay protocols may be adapted on other RT PCR platforms upon request

ELISA

TEST MENU

20

- BLOOD SCREENING
- COVID-19
- VIRAL HEPATITIS
- TORCH
- SEROLOGY
- EMERGING VIRUSES
- RETROVIRUSES
- AUTOIMMUNITY

20
22
23
25
26
28
30
31



BLOOD SCREENING

HEPATITIS B

CODE	KIT NAME	DESCRIPTION	TEST
SAG1ULTRA.CE.96	HBs Ag one version ULTRA	screening & mutants	96
SAG1ULTRA.CE	HBs Ag one version ULTRA	screening & mutants	192
SAG1ULTRA.CE.480	HBs Ag one version ULTRA	screening & mutants	480
SAG1ULTRA.CE.960	HBs Ag one version ULTRA	screening & mutants	960
BCAB.CE	HBc Ab	screening	96

HEPATITIS C

CODE	KIT NAME	DESCRIPTION	TEST
CVAB.CE.96	HCV Ab	version 4 - screening	96
CVAB.CE	HCV Ab	version 4 - screening	192
CVAB.CE.480	HCV Ab	version 4 - screening	480
CVAB.CE.960	HCV Ab	version 4 - screening	960

HUMAN IMMUNODEFICIENCY VIRUS (HIV) 1,2&0

CODE	KIT NAME	DESCRIPTION	TEST
IVCOMB.CE.96	HIV Ab&Ag	screening	96
IVCOMB.CE	HIV Ab&Ag	screening	192
IVCOMB.CE.480	HIV Ab&Ag	screening	480
IVCOMB.CE.960	HIV Ab&Ag	screening	960

HUMAN T-CELL LYMPHOTROPIC VIRUS (HTLV) I&II

CODE	KIT NAME	DESCRIPTION	TEST
HTLVABULTRA.CE.96	HTLV I&II Ab Version ULTRA	screening	96
HTLVABULTRA.CE	HTLV I&II Ab Version ULTRA	screening	192
HTLVABULTRA.CE.480	HTLV I&II Ab Version ULTRA	screening	480
HTLVABULTRA.CE.960	HTLV I&II Ab Version ULTRA	screening	960

BLOOD SCREENING



TREPONEMA PALLIDUM (SYPHILIS)

CODE	KIT NAME	DESCRIPTION	TEST
SIABULTRA.CE.96	Syphilis Ab Version ULTRA	screening	96
SIABULTRA.CE	Syphilis Ab Version ULTRA	screening	192
SIABULTRA.CE.480	Syphilis Ab Version ULTRA	screening	480
SIABULTRA.CE.960	Syphilis Ab Version ULTRA	screening	960
SIAB1ULTRA.CE.96	Syphilis Ab One Version ULTRA	screening	96
SIAB1ULTRA.CE	Syphilis Ab One Version ULTRA	screening	192
SIAB1ULTRA.CE.480	Syphilis Ab One Version ULTRA	screening	480
SIAB1ULTRA.CE.960	Syphilis Ab One Version ULTRA	screening	960

TRYPANOSOMA CRUZI (CHAGAS)

CODE	KIT NAME	DESCRIPTION	TEST
TCAB.CE.96	T.cruzi Ab	recAg screening	96
TCAB.CE	T.cruzi Ab	recAg screening	192
TCAB.CE.480	T.cruzi Ab	recAg screening	480
TCAB.CE.960	T.cruzi Ab	recAg screening	960

PLASMODIUM (MALARIA)

CODE	KIT NAME	DESCRIPTION	TEST
MALAB.CE.96	Malaria Ab	recAg screening	96
MALAB.CE	Malaria Ab	recAg screening	192
MALAB.CE.480	Malaria Ab	recAg screening	480
MALAB.CE.960	Malaria Ab	recAg screening	960

COVID-19

CODE	KIT NAME	DESCRIPTION	TEST
COV19M.CE	COVID-19 IgM	Determination of IgM antibodies	96
COV19M.CE.192	COVID-19 IgM	Determination of IgM antibodies	192
COV19A.CE	COVID-19 IgA	Determination of IgA antibodies	96
COV19A.CE.192	COVID-19 IgA	Determination of IgA antibodies	192
COV19G.CE	COVID-19 IgG	Determination of IgG antibodies	96
COV19G.CE.192	COVID-19 IgG	Determination of IgG antibodies	192
COV19GSPIKE.CE	COVID-19 Spike 1&2 IgG	Determination of IgG antibodies to COVID-19 Spike antigens	96
COV19GSPIKE.CE.192	COVID-19 Spike 1&2 IgG	Determination of IgG antibodies to COVID-19 Spike antigens	192
COV19CONF.CE	COVID-19 IgG Confirmation	Strips module Enzyme Immunoassay for the confirmation of IgG antibodies to COVID-19 major antigens	24
COV19TY.CE	COVID-19 IgG/IgM Confirmation and Typing	Strips module Enzyme Immunoassay for the confirmation of IgG and IgM antibodies to COVID-19 major antigens	24
ACE2-RBDNEUTR.CE	ACE2-RBD Neutralization Assay	Determination of inhibition activity of RBD-ACE2 binding induced by antibodies to SARS-CoV-2	96

VIRAL HEPATITIS



HEPATITIS A

CODE	KIT NAME	DESCRIPTION	TEST
AVAB.CE	HAV Ab	competitive	96
AVM.CE	HAV IgM	capture	96

HEPATITIS B

CODE	KIT NAME	DESCRIPTION	TEST
SAB.CE	HBs Ab	quantitative	96
SAGIULTRA.CE.96	HBs Ag one version ULTRA	screening & mutants	96
SAGIULTRA.CE	HBs Ag one version ULTRA	screening & mutants	192
SAGIULTRA.CE.480	HBs Ag one version ULTRA	screening & mutants	480
SAGIULTRA.CE.960	HBs Ag one version ULTRA	screening & mutants	960
SCONF.CE	HBs Ag Confirmation		20
SCONF.CE.40	HBs Ag Confirmation		40
BCAB.CE	HBc Ab	screening	96
BCM.CE	HBc IgM	capture quantitative	96
HBE.CE	HBe Ag&Ab		96

HEPATITIS C

CODE	KIT NAME	DESCRIPTION	TEST
CVAB.CE.96	HCV Ab	version 4 - screening	96
CVAB.CE	HCV Ab	version 4 - screening	192
CVAB.CE.480	HCV Ab	version 4 - screening	480
CVAB.CE.960	HCV Ab	version 4 - screening	960
CVM.CE	HCV IgM	sandwich	96
CCONF.CE	HCV Ab confirmation		12



VIRAL HEPATITIS

HEPATITIS D

CODE	KIT NAME	DESCRIPTION	TEST
DAB.CE	HDV Ab	total antibody	96
DAG.CE	HDV Ag		96
DIM.CE	HDV IgM	capture	96

HEPATITIS E

CODE	KIT NAME	DESCRIPTION	TEST
EVABULTRA.CE.96	HEV Ab Version ULTRA	screening	96
EVG.CE	HEV IgG	sandwhich	96
EVM.CE	HEV IgM	qualitative	96

TORCH



CYTOMEGALOVIRUS

CODE	KIT NAME	DESCRIPTION	TEST
CMVG.CE	CMV IgG	quantitative	96
CMVM.CE	CMV IgM	capture	96
CMVAV	CMV Antibody Avidity Test	NOT CE	48

TOXOPLASMA GONDII

CODE	KIT NAME	DESCRIPTION	TEST
TOXOG.CE	Toxo IgG	quantitative	96
TOXOM.CE	Toxo IgM	capture	96
TOXOAV	Toxo Antibody Avidity Test	NOT CE	48

RUBELLA VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
RUBG.CE	RUB IgG	quantitative	96
RUBM.CE	RUB IgM	capture	96
RUBAV	RUB Antibody Avidity Test	NOT CE	48

HERPES SIMPLEX VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
HSV1G.CE	HSV1 IgG	quantitative	96
HSV2G.CE	HSV2 IgG	quantitative	96
HSV1M.CE	HSV1 IgM	capture	96
HSV2M.CE	HSV2 IgM	capture	96
HSVG.CE	HSV1&2 IgG	quantitative	96
HSVM.CE	HSV1&2 IgM	capture	96



SEROLOGY

EPSTEIN BARR VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
VCM.CE	EBV VCA IgA		96
VCAG.CE	VCA IgG	quantitative	96
VCAM.CE	VCA IgM	capture	96
EBNG.CE	EBNA IgG	quantitative	96
EBNM.CE	EBNA IgM	quantitative	96
EAG.CE	Ea IgG		96
EAM.CE	Ea IgM		96
EBAV	EBV Avidity Test	NOT CE	48

HELICOBACTER PYLORI

CODE	KIT NAME	DESCRIPTION	TEST
HPG.CE	HP IgG	Pylori quantitative	96
HPM.CE	HP IgM	Pylori quantitative	96
HPA.CE	HP IgA	Pylori quantitative	96
CAGA.CE	CagA IgA	quantitative	96
CAGG.CE	CagA IgG	quantitative	96
HPAG.CE	HP Ag		48
K.ESTR	Extraction Kit for H. Pylori Ag		48

COXSACKIE B VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
COXBG.CE	CoxB IgG		96
COXBM.CE	CoxB IgM		96

SEROLOGY



MEASLES VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
MEAG.CE	Measles virus IgG		96
MEAM.CE	Measles virus IgM		96

MENINGITIS

CODE	KIT NAME	DESCRIPTION	TEST
MENG.CE	Meningitis IgG		96

YELLOW FEVER

CODE	KIT NAME	DESCRIPTION	TEST
YEVG.CE	Yellow Fever Virus IgG		96

TETANUS TOXOID

CODE	KIT NAME	DESCRIPTION	TEST
TETG.CE	TETOX IgG		96

TREPONEMA PALLIDUM (SYPHILIS)

CODE	KIT NAME	DESCRIPTION	TEST
SIM.CE	SYPH IgM	quantitative	96
SIABULTRA.CE.96	Syphilis Ab Version ULTRA	screening	96
SIABULTRA.CE	Syphilis Ab Version ULTRA	screening	192
SIABULTRA.CE.480	Syphilis Ab Version ULTRA	screening	480
SIABULTRA.CE.960	Syphilis Ab Version ULTRA	screening	960
SIAB1ULTRA.CE.96	Syphilis Ab One Version ULTRA	screening	96
SIAB1ULTRA.CE	Syphilis Ab One Version ULTRA	screening	192
SIAB1ULTRA.CE.480	Syphilis Ab One Version ULTRA	screening	480
SIAB1ULTRA.CE.960	Syphilis Ab One Version ULTRA	screening	960



SEROLOGY

CHLAMYDIA

CODE	KIT NAME	DESCRIPTION	TEST
CPA.CE	Chlamydia Pneumoniae IgA		96
CPG.CE	Chlamydia Pneumoniae IgG		96
CPM.CE	Chlamydia Pneumoniae IgM		96
CTA.CE	Chlamydia Trachomatis IgA		96
CTG.CE	Chlamydia Trachomatis IgG		96
CTM.CE	Chlamydia Trachomatis IgM		96

PARVO VIRUS B19

CODE	KIT NAME	DESCRIPTION	TEST
PARVOG.CE	Porvovirus B19 IgG		96
PARVOM.CE	Porvovirus B19 IgM		96

HUMAN PAPILLOMA VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
HPVG	Papilloma Virus IgG	R.U.O	96



EMERGING VIRUSES

TRYPANOSOMA CRUZI (CHAGAS)

CODE	KIT NAME	DESCRIPTION	TEST
TCAB.CE.96	T.cruzi Ab	recAg screening	96
TCAB.CE	T.cruzi Ab	recAg screening	192
TCAB.CE.480	T.cruzi Ab	recAg screening	480
TCAB.CE.960	T.cruzi Ab	recAg screening	960



EMERGING VIRUSES

PLASMODIUM (MALARIA)

CODE	KIT NAME	DESCRIPTION	TEST
MALAB.CE.96	Malaria Ab	recAg screening	96
MALAB.CE	Malaria Ab	recAg screening	192
MALAB.CE.480	Malaria Ab	recAg screening	480
MALAB.CE.960	Malaria Ab	recAg screening	960
MALTYPE.RUO	Malaria Ab Typing	R.U.O	48

WEST NILE VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
WNAB.CE	West Nile Virus Antibodies	screening	96
WNAB.CE.192	West Nile Virus Antibodies	screening	192
WNG.CE	West Nile Virus IgG		96
WNM.CE	West Nile Virus IgM		96

DENGUE VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
DENG.CE	Dengue virus IgG		96
DENM.CE	Dengue virus IgM		96
DENS1AG.CE	Dengue virus NS I Ag		96

ZIKA VIRUS

CODE	KIT NAME	DESCRIPTION	TEST
ZIKVM.CE	ZIKV IgM		96
ZIKVG.CE	ZIKV IgG		96
ZIKVAV.CE	ZIKV IgG Avidity Test		48

CHIKUNGUNYA

CODE	KIT NAME	DESCRIPTION	TEST
CHIKVG.CE	Chikungunya Virus IgG		96
CHIKVM.CE	Chikungunya Virus IgM		96



RETROVIRUSES

HUMAN IMMUNODEFICIENCY VIRUS (HIV) 1,2&0

CODE	KIT NAME	DESCRIPTION	TEST
IVCOMB.CE.96	HIV Ab&Ag	screening	96
IVCOMB.CE	HIV Ab&Ag	screening	192
IVCOMB.CE.480	HIV Ab&Ag	screening	480
IVCOMB.CE.960	HIV Ab&Ag	screening	960
IVP24	HIV p24 Antigen	NOT CE	96


HUMAN IMMUNODEFICIENCY VIRUS (HTLV) I&II

CODE	KIT NAME	DESCRIPTION	TEST
HTLVABULTRA.CE.96	HTLV I&II Ab Version ULTRA	screening	96
HTLVABULTRA.CE	HTLV I&II Ab Version ULTRA	screening	192
HTLVABULTRA.CE.480	HTLV I&II Ab Version ULTRA	screening	480
HTLVABULTRA.CE.960	HTLV I&II Ab Version ULTRA	screening	960
HTLVTYPE	HTLV I&II Ab typing	NOT CE	24

AUTOIMMUNITY



CODE	KIT NAME	DESCRIPTION	TEST
SSA52.CE	IgG anti SSA 52kD	rec Ag Quantitative	96
SSA60.CE	IgG anti SSA 60kD	rec Ag Quantitative	96
SSB.CE	IgG anti SSB	rec Ag Quantitative	96
RNP.CE	IgG anti U1 -sn RNP68	rec Ag Quantitative	96
DSDNA.CE	IgG anti dsDNA	plasmidic Ag Quantitative	96
SM.CE	IgG anti Sm	native Ag Quantitative	96
CENPB.CE	IgG anti Centromere B	rec Ag Quantitative	96
JO1.CE	IgG anti Jo-1	rec Ag Quantitative	96
SCL70.CE	IgG anti Scl-70	rec Ag Quantitative	96
ANAS.CE	ANA Screening IgG		96
ENAS.CE	ENA Screening IgG		96
ANA8PRO.CE	ANA 8 parameters profile		12
ENA6PRO.CE	ENA 6 parameters profile		12



Dia.Pro Diagnostic Bioprobes S.r.l.

Via G. Carducci, 27 • 20099 Sesto San Giovanni (Mi) • Italy
TEL + 39 02 27.00.71.61 /6450 • FAX + 39 02 44.38.67.71

@ info@diapro.it

P.IVA/C.F.: ITI I 924660159 - REG.IMP. I I 924660159
REA 1509959 - Capitale Sociale euro 50.000 i.v.